

Emerging Practice Summary: Reprocessing of N95 masks during severe shortages

Joe Cummings, PhD, Director, Vizient Technology Assessment Program

April 17, 2020

The CDC has recommended [strategies](#) for hospitals to optimize their supply of N95 masks stratified by the severity of their shortage. Contingency and crisis capacity strategies include options for [extended use and limited reuse](#) as well as use of other types of respirators, expired respirators and other controls.

The CDC has also provided [guidance](#) on the decontamination and reuse of respirators as a crisis strategy for hospital's experiencing severe shortages. Mask rotation using a minimum of 5 days under proper storage conditions is recommended as a passive decontamination strategy. Other promising methods noted by the CDC included ultraviolet germicidal irradiation (UVGI), hydrogen peroxide vapor (HPV) and moist heat. It is suggested that hospitals may want to primarily focus current efforts on these technologies because they have the most experience to-date. Many other decontamination technologies, however, have also been [proposed](#) and trialed.

Vizient has collected [polling data](#) from hospitals regarding N95 mask strategies. As of April 15, 78 hospitals reported reprocessing of N95 respirators. Forty-two respondents (54%) indicated they were using UVGI technology, 19% were using HPV, 8% were using heat, 6% were using ethylene oxide and 13% reported other unspecified methods. Reprocessing is being performed both in-house and through third-party reprocessors.

Ideally, the selection of disinfection technology should be provided by the respirator manufacturer after validation on their specific products, but third-party guidance may also be used in the absence of manufacturer recommendations according to CDC. At least one [manufacturer](#) of respirators has begun testing and approving various disinfection methods.

Important considerations for mask decontamination technologies include:

- efficacy of viral load deactivation
- maintenance of filtration efficiency
- no compromise of fit or weakening of any mask component

At this time, there are limitations in the evidence for each method. Most have not been tested directly against the SARS-CoV-2 virus; thus, they often extrapolate conclusions from use against different pathogens. Studies have also been conducted on only a limited number of mask types and using limited sample numbers. Since there are many versions of N95 masks, with different materials, layers, thicknesses and shapes, one method may work well for one mask type and not for another. Filtration efficiency will be affected by material degradation, but is also a function of induced changes in the electrostatic charge of the polypropylene strands. The number of decontamination cycles will also undoubtedly affect efficacy outcomes.

Further description of the pros and cons of the various mask decontamination technologies can be found at the N95decon consortium [website](#) and in the CDC [guidance](#). A bibliography is included at

the end of this document listing much of the primary evidence from published clinical studies derived from a search of the PubMed database. MedRxiv is another resource for some of the most recent, albeit non peer-reviewed, clinical evidence.

Findings from the reviewed evidence suggest UVGI has been studied on masks against various pathogens for more than 10 years and there is a fair amount of evidence in the clinical literature for its safety, efficacy and acceptance under controlled conditions. University of Nebraska Medical Center has published extensively on their technique using UVGI. Applied Research Associates has also published an FDA-sponsored study of UVGI for use on masks during public health emergencies that was conducted prior to the current pandemic.

Similarly, HPV had also undergone previous study by the Battelle Institute for use on respirators during public health emergencies. Duke University has also published their guidance on HPV disinfection of N95 respirators. The published research suggests HPV can be safely used specifically on N95 respirators with good viral deactivation and little compromise of materials over many cycles. In addition, the FDA has reviewed and granted emergency use authorization (EUA) for 4 different HPV-based technologies.

The number of cycles that a mask can be reprocessed is not well-established and will be highly dependent on the specific decontamination technology, operating parameters, mask model and findings from the quality assurance process conducted after every cycle. In general, for some specific regimens of UVGI and HPV there is data in the literature that support reprocessing up to 20 times without compromising filtration efficiency, straps or fit. FDA EUA labeling for select HPV technologies suggests a maximum between 2 and 20 cycles and anecdotal reports from some hospitals suggest they do not expect to reprocess masks more than about 5 to 7 times.

Conclusions

Reprocessing is a crisis strategy reserved for severe shortages when preferable alternatives are not available.

A number of hospitals have reported reprocessing strategies for masks due to crisis conditions.

The optimal reprocessing strategy has not been established due to limitations in the evidence.

There are a variety of options available that hospitals should consider based on their unique situations.

The selected decontamination method should effectively reduce the pathogen burden, maintain function and not affect the fit of the respirator.

Organizations should validate their process and outcomes to ensure safety.

Close inspection should be made after each decontamination cycle and affected masks should be discarded.

The number of reprocessing cycles should be minimized.

Due to fast-moving developments during this pandemic, efforts should be made to review any new evidence and recommendations regarding reprocessing when they become available.

Bibliographic Resources

- 3M. Decontamination Methods for 3M N95 Respirators-Revision 4 (April 2020). <https://multimedia.3m.com/mws/media/1824869O/decontamination-methods-for-3m-n95-respirators-technical-bulletin.pdf>. Accessed April 16, 2020.
- Anderegg L, Meisenhelder C, Ngooi CO, Liao L, Xiao W, Chu S, et al. A Scalable Method of Applying Heat and Humidity for Decontamination of N95 Respirators during the COVID-19 Crisis. *medRxiv* 2020.04.09.20059758; doi: <https://doi.org/10.1101/2020.04.09.20059758>.
- Baluja A, Arines J, Vilanova R, Bao-Varela C, Flores-Arias MT. UV Light Dosage Distribution over Irregular Respirator Surfaces. Methods and Implications for Safety. *medRxiv* 2020.04.07.20057224; doi: <https://doi.org/10.1101/2020.04.07.20057224>.
- Battelle Memorial Institute. Investigating Decontamination and Reuse of Respirators in Public Health Emergencies. FDA website. www.fda.gov/emergency-preparedness-and-response/mcm-regulatory-science/investigating-decontamination-and-reuse-respirators-public-health-emergencies. Accessed April 6, 2020.
- Battelle Memorial Institute. Final Report for the Bioquell Hydrogen Peroxide Vapor Decontamination for Reuse of N95 Respirators. www.fda.gov/media/136386/download. Accessed April 6, 2020.
- Bauchner H, Fontanarosa PB, Livingston EH. Conserving Supply of Personal Protective Equipment - A Call for Ideas. *JAMA*. 2020. <https://jamanetwork.com/journals/jama/fullarticle/2763590>.
- Bergman M, Viscusi DJ, Palmiero AJ, Powell JB, Shaffer RE. Impact of Three Cycles of Decontamination Treatments on Filtering Facepiece Respirator Fit. *J Inter Soc Resp Protection*. 2011;28(1):48-59.
- Bergman MS, Viscusi DJ, Heimbuch BK, Wander JD, Sambol AR, et al. Evaluation of Multiple (3-Cycle) Decontamination Processing for Filtering Facepiece Respirators. *J Engineered Fibers and Fabrics*. 2010;5(4):33-41.
- CDC. Recommended Guidance for Extended Use and Limited Reuse of N95 FFRs in Healthcare Settings. www.cdc.gov/niosh/topics/hcwcontrols/recommendedguidanceextuse.html. Accessed April 6, 2020.
- CDC. Decontamination and Reuse of Filtering Facepiece Respirators. www.cdc.gov/coronavirus/2019-ncov/hcp/ppe-strategy/decontamination-reuse-respirators.html. Accessed April 6, 2020.
- Card KJ, Crozier D, Dhawan A, Dinh M, Dolson E, Farrokhian N, et al. UV Sterilization of Personal Protective Equipment with Idle Laboratory Biosafety Cabinets During the Covid-19 Pandemic. *medRxiv* 2020.03.25.20043489; doi: <https://doi.org/10.1101/2020.03.25.20043489>.
- Cramer A, Tian E, Yu SH, Galanek M, Lamere E, Li J, et al. Disposable N95 Masks Pass Qualitative Fit-Test But Have Decreased Filtration Efficiency after Cobalt-60 Gamma Irradiation. *medRxiv* 2020.03.28.20043471; doi: <https://doi.org/10.1101/2020.03.28.20043471>.
- Derraik JGB, Anderson WA, Connelly EA, Anderson YC. Rapid evidence summary on SARS-CoV-2 survivorship and disinfection, and a reusable PPE protocol using a double-hit process. *medRxiv* 2020.04.02.20051409; doi: <https://doi.org/10.1101/2020.04.02.20051409>.
- Duan SM, Zhao XS, Wen RF, Huang JJ, Pi GH, et al. Stability of SARS coronavirus in human specimens and environment and its sensitivity to heating and UV irradiation. *Biomed Environ Sci*. 2003;16(3):246-255.
- Duke University. Decontamination and Reuse of N95 Respirators with Hydrogen Peroxide Vapor to Address Worldwide Personal Protective Equipment Shortages during the SARS-CoV-2 (COVID-19)

Pandemic. www.safety.duke.edu/sites/default/files/N-95_VHP-Decon-Re-Use.pdf. Accessed April 6, 2020.

FDA Guidance for Industry and FDA Staff: Enforcement Policy for Face Masks and Respirators During the COVID-19 Public Health Emergency (Revised April 2020). www.fda.gov/media/136449/download. Accessed April 6, 2020.

Fischer R, Morris DH, van Doremalen N, Sarchette S, Matson J, Bushmaker T, et al. Assessment of N95 respirator decontamination and re-use for SARS-CoV-2. *medRxiv* 2020.04.11.20062018; doi: <https://doi.org/10.1101/2020.04.11.20062018>.

Fisher EM, Shaffer RE. A method to determine the available UV-C dose for the decontamination of filtering facepiece respirators. *J Appl Microbiol*. 2011;110(1):287-295.

Fisher EM, Williams JL, Shaffer RE. Evaluation of microwave steam bags for the decontamination of filtering facepiece respirators. *PLoS One*. 2011;6(4):e18585.

Heimbuch BK, Harnish D. Research to Mitigate a Shortage of Respiratory Protection Devices During Public Health Emergencies -2019. Applied Research Associates website. www.ara.com/news/ara-research-mitigate-shortage-respiratory-protection-devices-during-public-health-emergencies. Accessed April 7, 2020.

Heimbuch BK, Wallace WH, Kinney K, Lumley AE, Wu CY, et al. A pandemic influenza preparedness study: use of energetic methods to decontaminate filtering facepiece respirators contaminated with H1N1 aerosols and droplets. *Am J Infect Control*. 2011;39(1):e1-9.

Lindsley WG, Martin SB, Jr., Thewlis RE, Sarkisian K, Nwoko JO, et al. Effects of Ultraviolet Germicidal Irradiation (UVGI) on N95 Respirator Filtration Performance and Structural Integrity. *J Occup Environ Hyg*. 2015;12(8):509-517.

Lore MB, Heimbuch BK, Brown TL, Wander JD, Hinrichs SH. Effectiveness of three decontamination treatments against influenza virus applied to filtering facepiece respirators. *Ann Occup Hyg*. 2012;56(1):92-101.

Kenney P, Chan BK, Kortright K, Cintron M, Havill N, Russi M, et al. Hydrogen Peroxide Vapor sterilization of N95 respirators for reuse. *medRxiv* 2020.03.24.20041087; doi: <https://doi.org/10.1101/2020.03.24.20041087>.

Kumar A, Kasloff S, Leung A, Cutts T, Strong JE, Hills K, et al. N95 Mask Decontamination using Standard Hospital Sterilization Technologies. *medRxiv* 2020.04.05.20049346; doi: <https://doi.org/10.1101/2020.04.05.20049346>.

Liao L, Xiao W, Zhao M, Yu X, Wang H, Wang Q, et al. Can N95 respirators be reused after disinfection? And for how many times? *medRxiv* 2020.04.01.20050443; doi: <https://doi.org/10.1101/2020.04.01.20050443>.

Mills D, Harnish DA, Lawrence C, Sandoval-Powers M, Heimbuch BK. Ultraviolet germicidal irradiation of influenza-contaminated N95 filtering facepiece respirators. *Am J Infect Control*. 2018;46(7):e49-e55.

N95Decon. A scientific consortium for data-driven study of N95 filtering facepiece respirator decontamination. www.n95decon.org/. Accessed April 6, 2020.

N95Decon Research Document (Version 1.2; 4/2/2020). Technical Report for Heat-Humidity-Based N95 Reuse Risk Management. https://static1.squarespace.com/static/5e8126f89327941b9453eeef/t/5e86d4b6c806bd1197556052/1585894582592/200402_N95DECON_Heat_technicalreport_v1.2_final.pdf. Accessed April 6, 2020.

N95Decon Research Document (Version 1.2; 4/2/2020). Technical Report for UV-C-Based N95 Reuse Risk Management.

https://static1.squarespace.com/static/5e8126f89327941b9453eeef/t/5e8541760211467623b1e4e4/1585791351199/200401_N95DECON_UV_technicalreport_v1.2_final.pdf. Accessed April 6, 2020.

N95Decon Research Document (Version 1.2; 4/2/2020). Technical Report for H2O2-Based N95 Reuse Risk Management.

https://static1.squarespace.com/static/5e8126f89327941b9453eeef/t/5e86d8be51cbbf1de3b3688e/1585895615567/200402_N95DECON_HP_V_technicalreport_v1.2_final.pdf. Accessed April 6, 2020.

Nazeeri AI, Hilburn I, Wu DA, Mohammed KA, Badal DY, Chan MHW, Kirschvink JL. An Efficient Ethanol-Vacuum Method for the Decontamination and Restoration of Polypropylene Microfiber Medical Masks & Respirators. *medRxiv* 2020.04.12.20059709; doi: <https://doi.org/10.1101/2020.04.12.20059709>.

Nebraska Medicine. N95 Filtering Facepiece Respirator Ultraviolet Germicidal Irradiation (UVGI) Process for Decontamination and Reuse. www.nebraskamed.com/sites/default/files/documents/covid-19/n-95-decon-process.pdf. Accessed April 6, 2020.

Nemeth C, Laufersweiler D, Polander E, Orvis C, Harnish D, et al. Preparing for an Influenza Pandemic: Hospital Acceptance Study of Filtering Facepiece Respirator Decontamination Using Ultraviolet Germicidal Irradiation. *J Patient Saf.* 2020.

Nogee D, Tomassoni A. Concise Communication: Covid-19 and the N95 Respirator Shortage: Closing the Gap [published online ahead of print, 2020 Apr 13]. *Infect Control Hosp Epidemiol.* 2020;1–4. doi: <https://doi.org/10.1017/ice.2020.124>.

Oral E, Wannomae KK, Connolly RL, Gardecki JA, Leung HM, Muratoglu OK, et al. Vapor H2O2 sterilization as a decontamination method for the reuse of N95 respirators in the COVID-19 emergency. *medRxiv* 2020.04.11.20062026; doi: <https://doi.org/10.1101/2020.04.11.20062026>.

Perkins DJ, Villescas S, Wu TH, Muller TB, Bradfute S, Foo-Hurwitz I, et al. COVID-19 Global Pandemic Planning: Decontamination and Reuse Processes for N95 Respirators. *medRxiv* 2020.04.09.20060129; doi: <https://doi.org/10.1101/2020.04.09.20060129>.

Salter WB, Kinney K, Wallace WH, Lumley AE, Heimbuch BK, et al. Analysis of residual chemicals on filtering facepiece respirators after decontamination. *J Occup Environ Hyg.* 2010;7(8):437-445.

Smith JS, Hanseler H, Welle J, Rattray R, Campbell M, Brotherton T, et al. Effect of various decontamination procedures on disposable N95 mask integrity and SARS-CoV-2 infectivity. *medRxiv* 2020.04.11.20062331; doi: <https://doi.org/10.1101/2020.04.11.20062331>.


Stanford Medicine COVID-19 Evidence Service. Addressing COVID-19 Face Mask Shortages [v1.3] Updated March 25, 2020. <https://stanfordmedicine.app.box.com/v/covid19-PPE-1-2>. Accessed April 6, 2020.

Viscusi DJ, Bergman MS, Eimer BC, Shaffer RE. Evaluation of five decontamination methods for filtering facepiece respirators. *Ann Occup Hyg.* 2009;53(8):815-827.

Viscusi DJ, Bergman MS, Novak DA, Faulkner KA, Palmiero A, et al. Impact of three biological decontamination methods on filtering facepiece respirator fit, odor, comfort, and donning ease. *J Occup Environ Hyg.* 2011;8(7):426-436.

Viscusi DJ, King WP, Shaffer RE. Effect of decontamination on the filtration efficiency of two filtering facepiece respirator models. *J Inter Soc for Resp Protection.* 2007;24:93-107.

 For further information, contact
Joe Cummings at
joe.cummings@vizientinc.com



DISCLAIMER: VIZIENT IS COMPILING INFORMATION AND EMERGING PRACTICES FROM MEMBERS TO AID IN KNOWLEDGE TRANSFER DURING THE COVID-19 RESPONSE. DECISIONS REGARDING WHETHER AND HOW TO UTILIZE ANY OF THESE PRACTICES SHOULD BE MADE BY HEALTH CARE PROVIDERS, AT THEIR OWN RISK, WITH CONSIDERATION OF INDIVIDUAL CIRCUMSTANCES. AS INFORMATION IS CHANGING RAPIDLY, VIZIENT ENCOURAGES YOU TO ALWAYS REFER TO THE CDC, YOUR STATE'S DEPARTMENT OF HEALTH, AND YOUR LOCAL PUBLIC HEALTH AUTHORITY FOR GUIDANCE. VIZIENT DOES NOT PROVIDE LEGAL, REGULATORY, OR MEDICAL ADVICE AND DISCLAIMS LIABILITY OR RESPONSIBILITY FOR THE ACCURACY, COMPLETENESS, AND/OR CLINICAL EFFICACY AND SAFETY FOR THE PRODUCTS OR PROCESSES CONTAINED HEREIN. MEMBERS SHOULD SEEK THEIR LEGAL COUNSEL'S ADVICE ON LOCAL, STATE, AND FEDERAL LEGAL/REGULATORY MATTERS. THE LINKS TO INFORMATION REFERENCED IN THIS DOCUMENT ARE THE PRODUCTS OF THE NAMED ORGANIZATIONS AND THEY ARE SOLELY RESPONSIBLE FOR THEIR CONTENT. FOR THE MOST UP-TO-DATE INFORMATION, PLEASE VISIT VIZIENT'S [DISASTER PREPAREDNESS PAGE](#). TO SUBMIT PRACTICES YOUR ORGANIZATION IS USING TO PREPARE FOR COVID-19, PLEASE EMAIL DISASTERRESPONSE@VIZIENTINC.COM.